

MeV Gamma Imaging by Fully Reconstructing Compton Scattering

A. Takada, T. Tanimori, H. Kubo, K. Miuchi, J. D. Parker, T. Mizumoto, Y. Mizumura, T. Sawano, Y. Matsuoka, S. Komura, S. Nakamura, M. Oda, S. Iwaki, K. Nakamura, S. Sonoda, D. Tomono (Kyoto Univ.)

- 1. Motivation & Detector concept
- 2. Performance
- 3. Applications
- 4. Summary



Detection of MeV gamma ray

Dominant process in MeV region -> Compton scattering

event

circle

gamma

- Elastic scattering between photon and electron.
- If detect momenta of scattered gamma ray and recoil electron
 We can obtain original direction and energy.

Compton Imaging



1st : interaction point & energy of recoil electron
 2nd : absorption point & energy of scattered gamma
 Not detect recoil direction

-> incomplete reconstruction

obtain source position by overlaying event circles artifacts appear in image



COMPTEL SMILE

Ryan, J. M., NewAR, 48, 199 (2004).

Electron-Tracking Compton Camera (ETCC)







A. Takada+, ApJ, 733(2011), 13

1st balloon experiment (SMILE-I)

Sub-MeV gamma-ray imaging Loaded-on-balloon Experiment

Launched on Sep. 1, 2006 @ Sanriku (ISAS/JAXA)

- Test flight using (10 cm)³ ETCC
- Measure diffuse cosmic and atmospheric gamma ray
 - 0.1 1 MeV, @ 35 km, 3 hours



ETCC

Measured : 420 events

Simulation : ~400 events (cosmic + atmospheric)

Compton kinematic test and Particle identify provided low-background observation.





ETCC for 2nd experiment

Target: Crab nebula

 5σ detection (40 km, several hours)

Requirements



Effective area $: > 0.5 \text{ cm}^2 (300 \text{ keV})$ Angular resolution : < 10° (600 keV) Sensitivity

 $: \times 100 \text{ SMILE-I}$

Improvements for SMILE-II

30 cm cube tracker

×√10 ×√10

- Updating of data acquisition system
- Improvement of imaging ability $\times 10$ Sensitivity will reach to (\times 100 SMILE-I)!







Event reconstruction



Back projection images





Angular resolution



Experiment 1: Confirmation of background rejection power



Experiment 1: Confirmation of background rejection power



Experiment 2: Observation of a week source

²²Na (40 kBq)

~2 M

Can ETCC detect gamma-ray source with low S/N?

Crab nebula : BG-gamma ≈ 0.01 : 1 Lead (2 mm)

Weak ²²Na -> ~100 ph/s come into ETCC 511 keV : BG = 0.02 : 1

Gamma-ray image has a clear excess. Significance of excess @ 511 keV is about 11σ during 5.5 h.



Experiment 3: Observation of time variation



Applications of ETCC

Medical imaging

- Tomography use as like PET/SPECT wide energy range

 -> possibility of new tracers
- Proton therapy monitoring of proton end point clear peak at high energy





Application of ETCC



Provide new efficient gamma-ray imaging detector for ¹³⁷Cs, ¹³⁴Cs in contaminated soils in Fukushima.

Project : Horiba & Kyoto-U & Canon supported by JST (Japan Science and Technology Agency)

- higher sensitivity
- compact and portable system
- smart data processing system and visualization software



10x10x15 cm³ ETCC with pressurized gas

At the first step

- soil bags were measured in Shirakawa city in Sep 2013.
 - Imaging contrasts
 - Separate spectrum components using images
 - Naïve estimation of radiation from energy spectra

very low dose < 0.05 μ Sv/h @ 1m (standard)

(0.6 ~ 2.0 $\mu\text{Sv/h}$ @ surface)



ETCC can measure absolute radiation dose.

Concept of environment monitoring was confirmed.

Now, we are developing a new camera, and planning some tests.

Summary

- We are developing an Electron-Tracking Compton Camera using a gaseous tracker.
- SMILE-II ETCC:
 - Effective area : ~1 cm² (< 300 keV)
 - Angular resolution : 5.3° (662 keV)
- ETCC has redundancies of background rejection
 - complete reconstruction using electron track
 - particles identify using dE/dx
 - Compton kinematic test using angle α
- Confirmation experiments:
 - detected gamma-ray source in high radiation field
 - detected a low S/N source
 - 511 keV, S/N = 0.02, live time = 2.0×10^4 s -> 10.5σ
 - confirmed the time variation sensitivity
 - 835 keV, 1 MBq, 3 m -> detected with 10 min at least
- Applications
 - Medical : proton therapy, SPECT, PET
 - Environment monitoring : in Fukushima, around reactor

Thank you for your attention!